

Feature	Continuous Monitoring	Routine & Ambulatory	Clinical Review ¹	Persyst Mobile ²	Multi-Patient Monitor	Advanced Review	Persyst Research
EEG Voltage Analysis	x	x	x			x	x
EEG Frequency Analysis	x	x	x			x	x
EEG Voltage & Freq. Mapping	x	x	x			x	x
Create Custom Frequency Filters	x		x			x	x
Export data in .CSV format	x		x			x	x
HFO Overlay							x
Batch processing							x