

Feature

RV: EEG Review	PAT: Read patient demographics	OLEV: Online event detection
DTEV: Outputs detection events to OEM database	RSEV: Reads recording system events	RSMNT: Reads recording system montages
DV: Reads synchronized digital video		
Alpha-Trace digital EEG	RV, PAT, OLEV, RSEV	*.alp
AstroMed Twin (Support Instructions)	RV, PAT	Logfile.qqq
Bio-logic BrainAtlas	RV, PAT	*.dat
Bio-logic Ceegraph V1	RV, PAT, RSEV	E*.dat
Bio-logic Ceegraph IV, V, VI, Vision and XL	RV, DV, PAT, OLEV, RSEV, RSMNT	*.eeg
BIOMECH (via Persyst Layout)	RV, PAT, OLEV, RSEV, RSMNT	*.lay
BioSemi BDF	RV, PAT	.bdf
BioSignal	RV, PAT, OLEV, RSEV	*.bpf
Blackrock Microsystems	RV, PAT	*.ns*
BrainLab	RV, DV, PAT	*.sig
Cadwell EasyWindows	RV, PAT, OLEV, DTEV, RSEV, RSMNT	*.eas
Cadwell Easy III and Arc EEG (Support Instructions)	RV, PAT, OLEV, DTEV, RSEV, RSMNT	*-1.ez3, *.arc
Cleveland Medical Devices Crystal-EEG	RV, PAT, OLEV, RSEV, RSMNT	*.crl
Compumedics Limited	RV, DV, PAT, OLEV, DTEV, RSEV, RSMNT	*.raw
Compumedics ProFusion	RV, PAT, OLEV, RSEV	*.sdy
DeltaMed (Support Instructions)	RV, PAT	*_0000.eeg
Deymed TruScan	RV, PAT	*.dat
DigiTrace	RV, PAT	*.w18,*.r27,*.w27
EBNeuro GaINT (Support Instructions)	RV, PAT, OLEV, RSEV	*.gnt
European Data Format 90	RV, PAT	*.edf
Grass Technologies AS40/Comet	RV, DV, PAT	Logfile.qqq
Grass Technologies AURA/AURA24	RV, DV, PAT	Logfile.qqq
Grass Technologies CTE	RV, DV, PAT	Logfile.qqq
Grass Technologies EDF	RV, PAT	*.edf
Grass Technologies TUFF	RV, PAT	*.ref
Grass Technologies TWin (Support Instructions)	RV, DV, PAT	Logfile.qqq
Lexicor	RV, PAT, OLEV, RSEV, RSMNT	*.dat
Lifelines/Trackit/EDF	RV, PAT, RSEV	*.edf, *.bdf
MayoVax	RV	*.eeg.1
MEF	RV	*.maf
Megsis Focus	RV	*.foc
Micromed	RV, PAT, OLEV	patient.info
Moberg CNS Monitor	RV, PAT, RSEV	*.trc
NEC Synafit	RV, PAT	*.dat
Nervus	RV, PAT, RSEV	*.e
Network Concepts Uniquant	RV, PAT	*.sif
NeuroScan	RV	*.cnt
Neuro Data NDF	RV, PAT	*.eeg
NeuroVirtual	RV	*.eeg*.psg
Nicolet Alliance NT	RV, DV, PAT, OLEV, RSEV, RSMNT	*.eeg
Nicolet BMSI 3.0+	RV	*.eeg
Nicolet Bravo	RV, DV, PAT, OLEV, RSEV, RSMNT	*.eeg
Nicolet Endeavor 3.x	RV, DV, PAT, OLEV, RSEV, RSMNT	*.xew
Nicolet NicoletOne cEEG	RV, DV, PAT, OLEV, RSEV, RSMNT	*.one
Nicolet NicoletOne nEEG/LTM/Monitor	RV, PAT, RSEV	*.e
Nicolet Satellite	RV, PAT, RSEV	*.dat
Nicolet Satellite ES	RV, PAT, RSEV	*.ebm
Nicolet Voyager	RV, PAT, OLEV, RSEV, RSMNT	data
NihonKohden 1100/2100/9100	RV, DV, PAT, OLEV, DTEV, RSEV, RSMNT	*.pnt
OSG BrainRT	RV, PAT, OLEV, RSEV	*-hdr.sig
Oxford 9200	RV	*.*
Oxford/Teca/Medelec DG32	RV, PAT, RSEV, RSMNT	*.mod
Oxford/Teca/Medelec Valor (Nervus)	RV, PAT, RSEV	*.e
Oxford/Teca/Medelec Profile	RV, PAT, RSEV	*.eeg
Polymate	RV, PAT, OLEV, RSEV	*.dat
QSI 9X00	RV, PAT	*.ano
Rembrandt	RV, PAT, OLEV, RSEV	*.msr
RMSIndia EEG-64	RV	*.eeg
Spectrum	RV, RSEV	*.00?
Stellate Harmonie (Support Instructions)	RV	*.sig
Stellate Monitor/Rhythm	RV, PAT	*.eeg
Telefactor CTE64	RV	*.dat
Telediagnostic (via Persyst Layout)	RV, PAT, OLEV, RSEV, RSMNT	*.lay
Vanguard	RV, PAT	b*
WalterGraphtek/Dantek PLEEG	RV, PAT, OLEV	*.wg1
Xltek	RV, DV, PAT, OLEV, RSEV, RSMNT	*.erd
Xltek EPWorks	RV, RSEV	*.iom